

DRUG CONTROL - Multiparameter, human-based**Catalogue Number:**

KG 1667 Drug Control - Level I
KG 1668 Drug Control - Level II
KG 1669 Drug Control - Level III

INTENDED USE

For use in the quality control of clinical chemistry assays, particularly drug residue analysis.

CHARACTERISTICS

BGT Drug Controls are based on lyophilised human serum. Constituent concentrations are available at 3 levels.

VALUE ASSIGNMENT

Each batch of serum is distributed to approximately 150 laboratories and values are assigned by a consensus of results obtained by these laboratories. A control range, for individual parameters and for each parameter method, is provided for each batch of serum. The control range is equivalent to the assigned mean \pm 2 S.D.

PREPARATION

Serum must only be reconstituted using the following procedure:

1. Open the vial carefully, avoiding any loss of material.
2. Reconstitute by pipetting exactly 5 ml of distilled water at +20 to +25°C into the vial.
3. Replace the rubber stopper and leave to stand for 30 minutes out of bright light before use.
4. Swirl gently several times during the reconstitution period to ensure that the contents are completely dissolved.
5. Prior to use, mix the contents by inverting the vials the formation of foam should be avoided. Ensure that no lyophilised material remains unreconstituted.
6. The serum is then ready for use with either a manual test or with an automated instrument.

STABILITY AND STORAGE

Unreconstituted serum is stable up to the expiry date shown on the side of each individual bottle when stored at +2 to +8°C. Once reconstituted the components of the Drug Control Sera are stable for 4 weeks at +2 to +8°C in the absence of bacterial contamination.

PRECAUTIONS AND WARNINGS

This serum has been tested for the HIV (Human Immunodeficiency Virus) antibody, ABsAg and HCV and found to be non-reactive. However, as no method can offer complete assurance as to the absence of infectious agents, this material should be handled as though capable of transmitting infectious disease.

This product has been developed for **IN VITRO** diagnostic use only.

Drug-o-trol Level 3 (TDM CONTROL 3)

Art.-Nr.: KG1669 Ch.-B.: 518DC

Inhalt 20 x 5ml

Verw. Bis: 2012-07

Bereich						
Parameter	Einheit	Zielwert	Von	Bis	Methoden	
Amikacin	µmol/l	47.3	37.8	56.8	Enzyme Immunoassay	
	µg/ml	27.7	22.1	33.3		
	µmol/l	47.0	37.6	56.4	Polarisation Fluoroimmunoassay	
	µg/ml	27.5	22.0	33.0		
Caffeine	µmol/l	92.2	73.8	111	HPLC (Reverse Phase)	
	µg/ml	17.9	14.3	21.5		
Carbamazepine	µmol/l	71.7	57.4	86.0	Enzyme Immunoassay	
	µg/ml	16.9	13.6	20.2		
	µmol/l	64.4	51.5	77.3	HPLC (Reverse Phase)	
	µg/ml	15.2	12.2	18.2		
	µmol/l	68.6	54.9	82.3	Polarisation Fluoroimmunoassay	
	µg/ml	16.2	13.0	19.4		
	µmol/l	59.4	47.5	71.3	Vitros 250/500/700/950/5.1 FS	
	µg/ml	14.0	11.2	16.8		
	µmol/l	66.1	52.9	79.3	Turbidimetric	
	µg/ml	15.6	12.5	18.7		
Cyclosporine	µmol/l	83.3	66.6	100	Roche Integra	
	µg/ml	19.7	15.7	23.7		
	µmol/l	76.9	61.5	92.3	KIMS	
	µg/ml	18.2	14.5	21.9		
	Cyclosporine	nmol/l	483	386	580	Polarisation Fluoroimmunoassay
		ng/ml	581	464	698	
	Digoxin	nmol/L	4.54	3.63	5.45	Vitros
		ng/ml	3.55	2.84	4.26	
nmol/L		4.72	3.78	5.66	Chemiluminescence	
ng/ml		3.69	2.95	4.43		
nmol/L		4.67	3.74	5.60	Enzyme Immunoassay	
ng/ml		3.65	2.92	4.38		
nmol/L		4.39	3.51	5.27	Polarisation Fluoroimmunoassay	
ng/ml		3.43	2.74	4.12		
nmol/L		4.01	3.21	4.81	Beckman CX/LX/Image	
ng/ml		3.13	2.51	3.75		
Ethosuximide	nmol/L	4.94	3.95	5.93	KIMS	
	ng/ml	3.86	3.08	4.64		
	nmol/L	4.81	3.85	5.77	Turbidimetric	
	ng/ml	3.76	3.01	4.51		
	Ethosuximide	µmol/l	947	758	1136	Enzyme Immunoassay
		µg/ml	134	107	161	
		µmol/l	865	692	1038	HPLC (Reverse Phase)
		µg/ml	123	98.1	148	
Gentamicin	µmol/l	17.2	13.8	20.6	Enzyme Immunoassay	
	µg/ml	7.96	6.39	9.53		

Drug-o-trol Level 3 (TDM CONTROL 3)

Art.-Nr.: KG1669 Ch.-B.: 518DC

Inhalt 20 x 5ml Verw. Bis: 2012-07

Bereich						
Parameter	Einheit	Zielwert	Von	Bis	Methoden	
Gentamicin	µmol/l	17.0	13.6	20.4	Polarisation Fluoroimmunoassay	
	µg/ml	7.87	6.30	9.44		
	µmol/l	17.8	14.2	21.4	Turbidimetric	
	µg/ml	8.24	6.57	9.91		
Lithium	µmol/l	14.4	11.5	17.3	KIMS	
	µg/ml	6.67	5.32	8.02		
	Lithium	mmol/l	1.93	1.70	2.16	Flame photometry
		mg/dl	1.34	1.18	1.50	
mmol/l		1.81	1.59	2.03	Ion selective electrode	
mg/dl		1.26	1.10	1.42		
Lithium (Vitros)	mmol/l	1.84	1.62	2.06	Spectrophotometric	
	mg/dl	1.28	1.12	1.44		
	Lithium (Vitros)	mmol/l	2.38	1.90	2.86	Vitros
		mg/dl	1.65	1.32	1.98	
Methotrexate	µmol/l	8.61	6.89	10.3	Enzyme Immunoassay	
	µg/ml	3.91	3.13	4.69		
	µmol/l	9.18	7.34	11.0	Polarisation Fluoroimmunoassay	
µg/ml	4.17	3.34	5.00			
	Paracetamol	mmol/l	1.49	1.19	1.79	Vitros
mg/l		225	180	270		
mmol/l		1.43	1.14	1.72	Colorimetric	
mg/l		216	172	260		
mmol/l		1.48	1.18	1.78	Polarisation Fluoroimmunoassay	
mg/l		224	179	269		
mmol/l		1.39	1.11	1.67	Enzymatic	
mg/l		210	168	252		
mmol/l	1.49	1.19	1.79	Turbidimetric		
mg/l	225	180	270			
Phenobarbitone	µmol/l	214	171	257	Enzyme Immunoassay	
	µg/ml	49.6	39.7	59.5		
	µmol/l	211	169	253	Polarisation Fluoroimmunoassay	
	µg/ml	49.0	39.2	58.8		
	µmol/l	191	153	229	HPLC (Reverse Phase)	
	µg/ml	44.3	35.5	53.1		
	µmol/l	204	163	245	Vitros	
	µg/ml	47.3	37.8	56.8		
µmol/l	216	173	259	Turbidimetric		
µg/ml	50.1	40.1	60.1			
Phenytoin	µmol/l	86.1	68.9	103	Vitros	
	µg/ml	21.7	17.4	26.0		
	µmol/l	90.2	72.2	108	Enzyme Immunoassay	
	µg/ml	22.8	18.2	27.4		
	µmol/l	85.3	68.2	102	HPLC (Reverse Phase)	
	µg/ml	21.5	17.2	25.8		
µmol/l	94.3	75.4	113	Polarisation Fluoroimmunoassay		
µg/ml	23.8	19.0	28.6			

Drug-o-trol Level 3 (TDM CONTROL 3)

Art.-Nr.: KG1669 Ch.-B.: 518DC Inhalt 20 x 5ml Verw. Bis: 2012-07

Bereich					
Parameter	Einheit	Zielwert	Von	Bis	Methoden
Phenytoin	µmol/l	100	80.0	120	Turbidimetric
	µg/ml	25.3	20.2	30.4	
	µmol/l	93.1	74.5	112	KIMS
	µg/ml	23.5	18.8	28.2	
Primidone	µmol/l	62.1	49.7	74.5	Polarisation Fluoroimmunoassay
	µg/ml	13.6	10.8	16.4	
	µmol/l	56.0	44.8	67.2	HPLC (Reverse Phase)
	µg/ml	12.2	9.78	14.6	
Salicylic Acid	mmol/l	2.84	2.27	3.41	Vitros
	mg/dl	39.2	31.4	47.0	
	mmol/l	2.81	2.25	3.37	Colorimetric Trinder
	mg/dl	38.8	31.1	46.5	
	mmol/l	2.92	2.34	3.50	Polarisation Fluoroimmunoassay
	mg/dl	40.3	32.3	48.3	
Theophylline	mmol/l	2.81	2.25	3.37	Enzymatic
	mg/dl	38.8	31.1	46.5	
	µmol/l	201	161	241	Vitros
	µg/ml	36.2	29.0	43.4	
	µmol/l	176	141	211	Enzyme Immunoassay
	µg/ml	31.7	25.4	38.0	
Tobramycin	µmol/l	179	143	215	Polarisation Fluoroimmunoassay
	µg/ml	32.3	25.8	38.8	
	µmol/l	178	142	214	HPLC (Reverse Phase)
	µg/ml	32.1	25.6	38.6	
	µmol/l	179	143	215	Turbidimetric
	µg/ml	32.3	25.8	38.8	
Valproic Acid	µmol/l	19.0	15.2	22.8	Enzyme Immunoassay
	µg/ml	8.89	7.11	10.7	
	µmol/l	18.7	15.0	22.4	Polarisation Fluoroimmunoassay
	µg/ml	8.75	7.02	10.5	
Vancomycin	µmol/l	17.4	13.9	20.9	Turbidimetric
	µg/ml	8.14	6.51	9.77	
	µmol/l	1009	807	1211	Enzyme Immunoassay
	µg/ml	146	116	176	
Valproic Acid	µmol/l	1014	811	1217	Polarisation Fluoroimmunoassay
	µg/ml	146	117	175	
	µmol/l	983	786	1180	Chemiluminescence
	µg/ml	142	113	171	
	µmol/l	1010	808	1212	Turbidimetric
	µg/ml	146	117	175	
Vancomycin	µmol/l	20.5	16.4	24.6	Enzyme Immunoassay
	µg/ml	30.5	24.4	36.6	
	µmol/l	19.6	15.7	23.5	Polarisation Fluoroimmunoassay
	µg/ml	29.1	23.3	34.9	
Vancomycin	µmol/l	19.3	15.4	23.2	Turbidimetric
	µg/ml	28.7	22.9	34.5	